


Spooky²

DH Experimental Frequencies for GeneratorX Offline



4 August 2019

By The Spooky² Team

Table of contents

DH Offline GX Presets Introduction	3
Installing the DH Offline GX Presets	4

DH Offline GX Presets Introduction

The *DH Experimental Frequencies* presets in the *Miscellaneous* folder of Spooky²'s Presets tab are quite popular. However, these frequencies only work for online mode of the Spooky² generators. The programs inside the factory DH presets are sometimes too large to fit inside the 200 frequency count limit of the GeneratorX program storage.

To get around this limitation of 200 frequencies per program on the GeneratorX, for the larger programs inside this collection I shortened the *Repeat Frequency* from 23 to a smaller value so that the program might fit inside of the frequency count limit. Even with this modification, some of the DH presets have so many different frequencies that some of these presets require more than one preset.

In the Presets notes of each preset, read the line that starts with *Offline Note*. This tells how long you need to run the preset to insure hitting each frequency 23 times. For the case of presets that have more than 1 offline preset for a DH online preset, the *Offline Note* tells you how much time to run each preset of the group of presets, and how much time it takes to run the entire set.



Most of the DH presets fit in a single offline preset. In this case, you load the single preset into the GeneratorX storage, and run the program (using the *Loop* button) for the time written in the *Offline Note* found in the *Preset Notes* of the preset. For the larger DH presets, you will need to load each program associated with a large preset into a different GeneratorX program location.

Installing the DH Offline GX Presets


You must first download the DH_GX_Offline_v2.zip archive from the Dropbox onto your computer. Extract all the contents of this archive.

Bring up a File Explorer window. 1. Click on *This PC* in the sidebar. 2. Double click the *Local Disk (C:)* drive.


Double click the *Spooky2* folder.



Double click the *Presets Collection* folder.



Double click the *User* folder.



In your *Downloads* folder in another *File Explorer* window, double click the archive. Now drag the folder from (1) to (2).

